

PHARMACEUTICAL

Terumo Syringe with and without Needle

Product Information Card

TERUMO syringes composed of barrel, plunger and gasket are sterile, single use medical devices intended for

- the injection of fluids, immediately after filling, into the human body
- aspiration of fluids
- sampling of blood or other body fluids.

The syringe portfolio covers all major size and tip types present on the general health care market

The finished products provide the required biocompatibility as set out in the applicable standards.

They do not contain any component made of PVC or Latex.

Components & raw materials of syringes

Component/part	Material type
1. Barrel	Polypropylene
2. Plunger	Polystyrene (1ml-size) or Polypropylene
3. Gasket	Thermoplastic elastomer
Lubricant	Silicone Oil
Blister	Polymer film



Regulatory status

CE certified



Manufacturer	TERUMO (PHILIPPINES) CORPORATION (TPC) 124 East Main Avenue / Laguna Technopark, Binan, Philippines
Notified Body	TÜV Rheinland LGA Products GmbH, Germany 
EC Representative	Terumo Europe N.V. Interleuvenlaan 40 3001 Leuven, Belgium

Please always consult the product label for a complete overview of warnings, cautions and precautions prior to actual use. Do not use if the packaging is damaged. Please note that some products may not be available in your region, please reach out to your local Terumo Medical Care Solutions - Pharmaceutical contact for details. All brand names are trademarks or registered trademarks owned by TERUMO CORPORATION, its affiliates, or unrelated third parties.

Syringes without needle	Reference	Nominal Capacity / ml	Graduated Capacity / ml	Increments / ml	Tip Type	Manufacturer	Unit Box Quantity	Shipping Carton Quantity	Grossweight / g (Shipping Carton)
	SS+01T1	1	1	0.01	Luer-slip, centric	TPC	100	1800	7.355
	SS+03L1	3	3	0.1	Luer-lock, centric	TPC	100	1800	7.539
	SS+05L1U	5	6	0.2	Luer-lock, centric	TPC	100	1200	7.488
	SS+10ES1	10	12	0.2	Luer-slip, eccentric	TPC	100	1200	11.827
	SS+10L1	10	12	0.2	Luer-lock, centric	TPC	100	1200	11.875
	SS+20ES1	20	25	1.0	Luer-slip, eccentric	TPC	50	600	9.964
	SS+20L1	20	25	1.0	Luer-lock, centric	TPC	50	600	9.877
	SS+30L1	30	35	1.0	Luer-lock, centric	TPC	25	200	6.745
	SS+50L1	50	60	1.0	Luer-lock, centric	TPC	25	100	5.070
	SS+50ES1	50	60	1.0	Luer-slip, eccentric	TPC	25	100	5.086
Syringes with needle	Reference	Nominal Capacity / ml	Graduated Capacity / ml	Increments / ml	Needle Dimensions (Outer Diameter / Length)	Manufacturer	Unit Box Quantity	Shipping Carton Quantity	Grossweight / g (Shipping Carton)
	SS+01T25161	1	1	0.01	25G (0.5mm) / 16mm (5/8")	TPC	100	1800	9.603
	SS+01T26131	1	1	0.01	26G (0.45mm) / 13mm (1/2")	TPC	100	1800	9.600
	SS+05S21381	5	6	0.2	21G (0.8mm) / 38mm (1 1/2")	TPC	100	1200	9.262
	SS+05S23321	5	6	0.2	23G (0.6mm) / 32mm (1 1/4")	TPC	100	1200	9.249

Selected portfolio – other specifications may be available too. Lead time to be confirmed.

Please note that some items may not be available in all regions, check with your local Medical Care Solutions - Pharmaceutical contact for details.

All brand names are trademarks or registered trademarks owned by TERUMO CORPORATION, its affiliates, or unrelated third parties.



Terumo Europe N.V.
Phone: +32 16 38 12 11

Terumo Corporation
Phone: +81 3 3374 81 8111

Terumo Medical Corporation
Phone: +1 732 302 4900
Toll free: +1 800 283 7866

www.terumopharmaceuticalsolutions.com