


Quality time for better care



PHARMACEUTICAL

SurGuard™3 Safety Hypodermic Needle



Safety made easy

- Supports all 3 ways of activation
- Audible, visible and tangible confirmation²
- User comfort – leading to patient's confidence²
- EN ISO 23908 compliant

Drug Delivery Devices

Features and benefits

- Hinged safety sheath design – stable, reliable, easy to adopt²
- Sharper needles leading to decreased penetration resistance¹
- Sheath locks to both the needle and the hub to provide added security and visible, audible and tactile feedback²
- Finger and thumb grips designed to suit all hand sizes
- Can be used on luer slip tip and luer lock syringes
- 3 ways of shield activation to suit individual clinician needs and preferences: hard surface or finger or thumb

Sheath activation



Images are for visualization purpose only.

¹ Comparative Hypodermic Needle Sharpness Among Major Brands: Results of an Independent Study 2009.

² Data on file

Refer to the product labels and instructions for use for warnings, cautions and precautions.

Please contact us for more details.

SurGuard™3 Safety Hypodermic Needle

References	Hub Color	Diameter (OD)		Length	
		Gauge	mm	mm	Inch
SG3+1813	Pink	18	1.2	13	1/2
SG3+1825	Pink	18	1.2	25	1
SG3+1838	Pink	18	1.2	38	1 1/2
SG3+1916	Ivory	19	1.1	16	5/8
SG3+1938	Ivory	19	1.1	38	1 1/2
SG3+2025	Yellow	20	0.9	25	1
SG3+2038	Yellow	20	0.9	38	1 1/2
SG3+2051	Yellow	20	0.9	51	2
SG3+2125	Green	21	0.8	25	1
SG3+2138	Green	21	0.8	38	1 1/2
SG3+2151	Green	21	0.8	51	2
SG3+2213	Black	22	0.7	13	1/2
SG3+2225	Black	22	0.7	25	1
SG3+2238	Black	22	0.7	38	1 1/2
SG3+2325	Blue	23	0.6	25	1
SG3+2338	Blue	23	0.6	38	1 1/2
SG3+2516	Orange	25	0.5	16	5/8
SG3+2525	Orange	25	0.5	25	1
SG3+2609	Brown	26	0.45	9	3/8
SG3+2713	Grey	27	0.4	13	1/2
SG3+3013	Yellow	30	0.30	13	1/2

Please note that some products may not be available in all regions, please reach out to your local Terumo Medical Care Solutions - Pharmaceutical contact for details. All brand names are trademarks or registered trademarks owned by TERUMO CORPORATION, its affiliates, or unrelated third parties.

CE 0197

Sterilized using irradiation.
The device has no components made of natural rubber latex or PVC.

SurGuard™3
Safety Hypodermic Needle

TERUMO

Terumo Europe N.V.
+32 16 38 12 11

Terumo Corporation
+81 3 6742 8500

Contact details are subject to change.
Please consult our website.

Terumo Medical Corporation
Phone: +1 732 302 4900
Toll free (U.S.A.): +1 800 283 7866

www.terumopharmaceuticalsolutions.com