

# Terumo Syringe

*3-part Syringe (with or without Needle)*



## Product Solutions

Your **syringes of choice**

**NEW**



 **TERUMO**

**Helping you  
to help people**



Offering  
**sustainable  
solutions**

# Disposable syringes

## for single use

The syringes are designed to be as safe as possible, convenient and represent the current state of the art in modern injection technology.

Carefully selected materials – not made with PVC or natural rubber latex – together with a controlled production process result in products of excellent quality and performance.



Terumo Syringes and  
Terumo Syringes with Needles

The range includes

- Concentric Luer Slip tip
- Eccentric Luer Slip tip
- Luer Lock tip



# Smooth and perfectly controlled injection for reliable drug delivery

Made to support you in your daily practice

Every single detail is the result of Terumo's long time experience in this field which is appreciated by users throughout the world.

**Performance enhancing features:**

- siliconisation for optimal and smooth plunger mobility
- the integral gasket ensures a firm and distinct end-stop, this helps to prevent the plunger from accidentally being pulled out
- finger grips to support easy handling

**Safety enhancing features:**

- transparent barrel provides easy recognition of air bubbles
- clear graduated scale allows accurate dosage



**ISO colour-coding of the needles:**

- quick identification of the required size
- easy recognition of air bubbles or flashback through the translucent hub


# Choose from various sizes optionally with needles

## Terumo Syringes

Sterile 3-part syringes  
in a choice of frequently  
needed volumes

## Terumo Syringes with needles

Sterile 3-part syringes  
pre-connected with a sterile  
hypodermic needle



### Syringe in 1 ml volume

available with these needles:



**25G Cannula (ø 0.5 mm)**

Length 16 mm, Orange

**26G Cannula (ø 0.45 mm)**

Length 13 mm, Brown

**27G Cannula (ø 0.4 mm)**

Length 13 mm, Grey

### Syringe in 2.5/3 ml volume

available with these needles:



**21G Cannula (ø 0.8 mm)**

Length 16/25/38 mm, Green

**22G Cannula (ø 0.7mm)**

Length 38 mm, Black

**23G Cannula (ø 0.6 mm)**

Length 25/32 mm, Blue

**25G Cannula (ø 0.5 mm)**

Length 16 mm, Orange

### Syringe in 5 ml volume

available with these needles:



**21G Cannula (ø 0.8 mm)**

Length 38 mm, Green

**22G Cannula (ø 0.7mm)**

Length 38 mm, Black

**23G Cannula (ø 0.6 mm)**

Length 32 mm, Blue

### Syringe in 10 ml volume

available with these needles:



**20G Cannula (ø 0.9 mm)**

Length 38 mm, Yellow

**21G Cannula (ø 0.8 mm)**

Length 38 mm, Green

**22G Cannula (ø 0.7mm)**

Length 38 mm, Black

For further ordering information, please refer to the Terumo Syringes /  
Terumo Syringes with needle product information card.

CE  
0197

Does not contain  
components made of PVC  
or natural rubber latex.





## Dedicated to making things better

**Terumo Syringes and Terumo Syringes with Needles will enhance Terumo's advanced Injection & Access portfolio**

We understand our devices are an essential element for treating patients and saving lives – especially in crisis impacted scenarios. It is therefore our commitment to ensure stability of supply whilst bringing new products to the market.

With the new Terumo Syringes we provide solutions to medical institutions and caregivers.

## Contributing to Society through Healthcare



**Terumo Corporation**  
+81 3 6742 8500

**Terumo Europe N.V.**  
+32 16 38 12 11

Your Terumo partner