

Quality time for better care



PHARMACEUTICAL

SurGuard™ 2



The convenient way to prevent needle stick injuries




















- Easy-to-learn One-Hand-Activation¹
- Visual, audible and tactile confirmation¹
- Full advantage of Terumo's T-SHARP™ Needle Technology²
- EN ISO 23908 compliant

With the development of the SurGuard™2 shield and the full integration of Terumo's outstanding experience in cannulas for hypodermic injections, we offer a convenient solution to enable a safe work environment:

- Limited training demand thanks to the intuitive design and user experience-supported Instructions for Use.¹
- Applicable for all major injection techniques: subcutaneous, intramuscular and intravenous (at various angles).
- User confidence is well supported by the provision of visual, audible and tactile feedback on successful activation.¹
- Sheath position pre- and post-injection can be chosen according to the user's preference, allowing an unimpaired view of the puncture point and requiring little interaction to activate the safety feature.
- All the advantages of Terumo's T-SHARP™ needle technology further contribute to user and patient comfort: low penetration and almost linear drag forces (virtually no knocking effect), with the aim of reducing pain as much as possible.²
- Well-balanced range of diameters and lengths to allow selection of the most suitable specification for any drug-application combination.

Our long presence in numerous markets and the increasing demand in occupational health support have led to a broad global registration coverage – and our teams are ready to extend it for you!

SurGuard™2 Safety Hypodermic Needle

References	Hub Color		Diameter (OD)		Length	
			Gauge	mm	mm	Inch
SG2+1825	Pink		18	1.2	25	1
SG2+1838	Pink		18	1.2	38	1 1/2
SG2+1925	Ivory		19	1.1	25	1
SG2+1938	Ivory		19	1.1	38	1 1/2
SG2+2025	Yellow		20	0.9	25	1
SG2+2038	Yellow		20	0.9	38	1 1/2
SG2+2125	Green		21	0.8	25	1
SG2+2138	Green		21	0.8	38	1 1/2
SG2+2225	Black		22	0.7	25	1
SG2+2238	Black		22	0.7	38	1 1/2
SG2+2325	Blue		23	0.6	25	1
SG2+2338	Blue		23	0.6	38	1 1/2
SG2+2516	Orange		25	0.5	16	5/8
SG2+2525	Orange		25	0.5	25	1
SG2+2538	Orange		25	0.5	38	1 1/2
SG2+2609	Brown		26	0.45	9	3/8
SG2+2613	Brown		26	0.45	13	1/2
SG2+2713	Grey		27	0.4	13	1/2
SG2+3013	Yellow		30	0.3	13	1/2

Shipping carton: 800 pieces (8 x unit box)
Sterilisation with irradiation

CE 0197 

The device has no components made of PVC/DEHP.

The device has no components made of natural rubber latex.

Recommendations for use



Refer to the product labels and instructions for use for warnings, cautions and precautions.

¹ Data on file

² Comparative Hypodermic Needle Sharpness Among Major Brands, 2009

Please note that some products may not be available in all regions. Check with your local Terumo Medical Care Solutions - Pharmaceutical contact for details. All brand names are trademarks or registered trademarks owned by TERUMO CORPORATION, its affiliates, or unrelated third parties.



Terumo Europe N.V. +32 16 38 12 11
Terumo Corporation +81 3 6742 8500

Terumo Medical Corporation
Phone: +1 732 302 4900
Toll free (U.S.A.): +1 800 283 7866

Contact details are subject to change.
Please consult our website.

www.terumopharmaceuticalsolutions.com

