

Quality time for better care



PHARMACEUTICAL

PLAJEX™

Ready-to-Fill Polymer Syringe with Staked Needle



Designed to meet the needs
of a variety of drug types



PLAJEX™ is designed to meet the needs of a wide variety of drug types, including those requiring lower interaction with containers.

PLAJEX™ integrates a staked needle with a robust 1 ml long silicone oil-free polymer syringe. This solution comes nested and sterilized, ready for filling and finishing. Offering excellent strength and clarity, and Terumo's proprietary i-coating™ technology on the plunger stopper, the silicone oil-free syringe system is designed to preserve sensitive biotherapeutics¹

PLAJEX™ offers the quality and supply reliability needed to meet parenteral delivery challenges and product commercialization goals

- Ready-to-fill format
- Silicone oil-free
- Tungsten-free
- Adhesive-free
- Higher resistance to breakage
- Autoclave advantage

i-coating™

- Terumo's proprietary i-coating™ is applied to the surface of the stopper and bonded to the substrate. It is not a laminate and it remains flexible with the substrate.
- Provides a silicone oil-free container for drug products requiring lower interaction with containers.
- This thin surface coating provides consistent and predictable gliding force for the syringe producing a silicone oil-free solution.

Images are for visualization purpose only.

Please note that some items may not be available in all regions, please reach out to your local Terumo Medical Care Solutions - Pharmaceutical contact for details. All brand names are trademarks or registered trademarks owned by TERUMO CORPORATION, its affiliates, or unrelated third parties.

PLAJEX™
Ready-to-Fill Polymer Syringe



¹ Functional Evaluation and Characterization of a Newly Developed Silicone Oil-Free Prefillable Syringe System, Yoshino et al (2014) J Pharm Sc, 103: 1520-1528



Terumo Europe N.V.
+32 16 38 12 11

Terumo Corporation
+81 3 6742 8500

Contact details are subject to change.
Please consult our website.

Terumo Medical Corporation
Phone: +1 732 302 4900
Toll free (U.S.A.): +1 800 283 7866

www.terumopharmaceuticalsolutions.com